Description: This antiserum was raised in a guinea pig by immunization with a bovine serum albumin (BSA) conjugate of synthetic bombesin analogue, N^{α} -glycyl -[Gln¹]-bombesin¹). The product vial contains 50 μ L of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with bombesin.

Immunogen: Synthetic N^{α} -glycyl- $[Gln^{\dagger}]$ -bombesin-BSA conjugate **Host:** Guinea pig

Amino Acid Sequence of Bombesin²⁾: pEQRLGNQWAV GHLM-NH2

Product Form: Lyophilized unpurified serum Size: 50 μ L

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

Suggested Working Dilution Range: 1:1,000 (final dilution ~1:7,000) for radioimmunoassay; 1:200-1,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay): Bombesin 100%, N^{α} -glycyl-[Gln¹]-bombesin 100%, N^{α} -tyrosyl-[Gln¹]-bombesin 100%, des-14-amide-bombesin 25%, bombesin (4-14) 90%, bombesin (6-14) 90%, bombesin (8-14) 0%, ranatensin C 2.6%, GRP (porcine) < 0.001%¹⁾.

Positive Control (immunohistochemistry): Rat duodenum

Species Tested: Guinea pig, dog, porcine, monkey

REFERENCES:

- 1) C. Yanaihara, N. Yanaihara et al., Bombesin-like immunoreactivity in mammalian tissues. Biomedical Research 1: 96-100, 1980
- 2) A. Anastasi, V. Erspamer and M. Bucci, Isolation and structure of bombesin and alytesin of the European amphibians, Bombina and Alytes. Expreimentia, 27: 166-167, 1971

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Yanaihara Instituta Inc.

2480-1 Awakura, Fujinomiya, Shizuoka, 418-0011 Japan

TEL: +81-544-22-2771 FAX:+81-544-22-2770