

## **Anti PTH (1-15) (Mouse, Rat) Serum**

Cat. No. Y421

Lot No. 274141024

**Description:** This antiserum was raised in a rabbit by immunization with a keyhole limpet hemocyanin (KLH) conjugate of synthetic PTH (1-15) (mouse, rat) peptide. The product vial contains 50  $\mu$ L of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry, or any other immunoreaction with PTH (mouse, rat).

**Immunogen:** Synthetic PTH (1-15) (mouse, rat)-KLH conjugate

**Host:** Rabbit

**Amino Acid Sequence of PTH (1-15) (Mouse, Rat)<sup>1)</sup>:** 1 15  
AVSEIQLMHN LGKHL

**Product Form:** Lyophilized unpurified serum

**Size:** 50  $\mu$ L

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:3,000 (final dilution ~1:21,000) for radioimmunoassay; 1: 1,000-4,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on enzymeimmunoassay): PTH (1-15) (mouse, rat) 100%, PTH (1-34)-NH<sub>2</sub> (mouse, rat) 10%

**Positive Control** (immunohistochemistry): Rat parathyroid gland

**Species Tested:** Rat

### **REFERENCES:**

1) G. heinrich, H.M. Kronenberg et al., Gene encoding parathyroid hormone. Nucleotide sequence of the rat gene and deduced amino acid sequence of rat pre- parathyroid hormone. Journal of Biological Chemistry 259:3320, 1984

**FOR RESEARCH LABORATORY USE ONLY**  
DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

  
Yanaihora Institute Inc.

2480-1 Awakura, Fujinomiya, Shizuoka, 418-0011 Japan  
TEL: +81-544-22-2771 FAX: +81-544-22-2770

For more information, please visit our internet website: <http://www.yanaihora.co.jp> or contact us to: [ask@yanaihora.co.jp](mailto:ask@yanaihora.co.jp)